



FORD TOURNEO CONNECT 1.5 TDCI TITANIUM 120PS [PO4206]

£23,995.00 ~~£22,495.00~~ Excluding VAT



DETAILS

Headroom (in.)	59
Wheelchair Armrest Width (in.)	26
Transmission	<u>Automatic</u>
Mileage	32656
Seats	4 + Wheelchair
Model Year	2017
Colour	Deep Impact Blue
Engine	1.5 TDCi EcoBlue
Fuel Type	Diesel
Power (PS)	120

This used Ford Tourneo Connect WAV features our inclusive Brotherwood conversion with two rear seats, 59" internal headroom and central wheelchair location. Finished in popular Deep Impact Blue this popular small WAV also features our PAWRS powered restraint system and STORQ system for fast, easy loading.

Discount subject to delivery before 30th June 2026

SKU: PO4206

We're proud
to be an
RAC Approved
Dealer

RAC
Approved Dealer



VEHICLE DESCRIPTION

ULEZ Compliant

Vehicle Summary

Dimensions

	mm	inches
Entry Height	1427	56.2
Max. Internal Headroom	1500	59
Lowered Floor Length	1511	59.5
Lowered Floor Width	790	31.1
Access Ramp Length	1110	43.7
Access Ramp Width	735	28.9
Access Ramp Angle		12.3°
Exterior Height	1852	73
Exterior Length	4418	174
Exterior Width	2137	84

Every effort has been made to ensure the accuracy of all information. However, errors may occur. Please check with the sales team about items which may affect your decision to purchase. The majority of our customers are exempt VAT, hence all prices shown exclude VAT.



Conversion Features

- Four Point Webbing Restraint
- High-Traction Ramp Surface
- PAWRS Power Assisted Wheelchair Restraint System
- STORQ Self-Tensioning Restraints
- Privacy Glass
- Counterbalanced Access Ramp
- Lowered, Level Floor

Specification Includes

- Auto wipers and headlamps
- DAB Radio/CD with SYNC
- Drivers 8-way seat, with lumbar support
- Quickclear heated windscreen & washer jets
- Manual air conditioning
- Cruise Control with Adjustable Speed Limiter
- Front parking sensors
- Rear parking sensors
- Panorama roof with Privacy Glass



GET A PART EXCHANGE VALUATION

Enter your details below to get a valuation for your vehicle in part exchange.

ASK US A QUESTION

FREEPHONE:

[\(+44\) 0330 1747 555](tel:+4417475550330)

Contact us today for more information and to discuss your needs. Our dedicated, experienced team are on hand to answer all of your questions about Wheelchair Accessible Vehicles.