



# FORD TOURNEO CONNECT 2.0L ACTIVE ECOBLUE 122PS [S6521]

£40,847.51 ~~£36,847.51~~ Excluding VAT

## DETAILS



Headroom (in.)	56
Wheelchair Armrest Width (in.)	25.5
Transmission	Manual
Mileage	124
Seats	4 + Wheelchair
Model Year	2025
Colour	Stardust Silver
Engine	2.0 EcoBlue
Fuel Type	Diesel
Power (PS)	122

This brand-new Ford Tourneo Connect comes fully equipped with 'Active' line specification including navigation, Apple CarPlay / Android Auto, and front & rear parking sensors. Currently in final production, reserve this new model today to choose your seating layout.

Standard conversion specification includes our signature flat, level lowered floor, PAWRS combined winch & restraint system and STORQ self-tensioning rear restraint points, for effortless access and fast, safe loading.

**SKU:** S6521





# VEHICLE DESCRIPTION

ULEZ Compliant

## Vehicle Summary

### Dimensions

	mm	inches
Entry Height	1395	55
Max. Internal Headroom <sup>1</sup>	1420	56
Lowered Floor Length	1640	64.5
Lowered Floor Width at Entry	760	30
Lowered Floor Width at Seats	705	27.75
Max. Wheelchair Armrest Width <sup>2</sup>	645	25.5
Access Ramp Length	1130	44.5
Access Ramp Width	730	28.75
Access Ramp Angle		12.5°
Exterior Height	1832	72
Exterior Length	4500	177
Exterior Width (inc. Mirrors)	2100	82.5

<sup>1</sup> Internal headroom may vary depending on wheelchair travelling position.

<sup>2</sup> In standard 2 rear seat configuration.

Every effort has been made to ensure the accuracy of all information. However, errors may occur. Please check with the sales team about items which may affect your decision to purchase. The majority of our customers are exempt VAT, hence all prices shown exclude VAT.



#### Conversion Features

- PAWRS Power Assisted Wheelchair Restraint System
- 2 Rear Seats
- STORQ Self- Tensioning Restraint System
- Four Point Webbing Restraint
- High-Traction Ramp Surface
- Counterbalanced Access Ramp
- Lowered, Flat, Level Floor

#### Specification Includes

- Adaptive Cruise Control Pack: Lane Keeping Aid, Pre-Collision Assist, Driver Impairment Monitor, Front & Rear Parking Aids, Hill Start Assist, Stop & Go (Auto only) Blind Spot Assist with Cross Traffic Alert and Exit Warning, Rear View Camera



- Android Auto and Apple CarPlay with USB or Wireless Connectivity
- DAB Radio with Bluetooth, voice control and Navigation (includes 1 year free subscription to Connected Navigation features)
- Dual Power Fold Heated Mirrors
- Electric Parking Brake with Hill Start Assist
- Front Fog Lamps
- Halogen Headlamps
- Keyless Start with Double Locking
- Quickclear Electric Windshield Defroster
- Rear Dark Tinted Glass
- Voice Control
- 17" Premium Painted Alloy Wheel
- Acoustics Package High – Improved Firewall and Rear Wheelarch insulation
- Auto Headlamps On/Off
- Body Coloured Plastic Front Bumper & Skid Plates
- Driver Assistance Package: Pre-Collision Assist with Auto Emergency Braking (includes Pedestrian and Cyclist Detection), Forward Collision Warning, Evasive Steer Assist, Intersection Assist, Cruise Control with Adjustable Speed Limiter, Lane-Keeping System with Driver Alert, Front and Rear Parking Sensors
- Driver heated seat with lumbar support, under seat stowage and manual 4-way adjuster. Single passenger heated seat with lumbar support, under seat stowage and 4-way manual adjuster and Fold Flat Passenger Seat. Contains Driver and Passenger seat airbags and Curtain airbags
- Front Parking Aid
- ICE Feature Pack 4: DAB radio with 10" touchscreen with voice and gesture control and navigation, Apple CarPlay, Android Auto. 2 USB type-C ports, Bluetooth® and 6 speakers
- Keyless Start with Double Locking
- Manual Air Conditioner
- Reverse Parking Aid
- Wheelarch Mouldings – Black

---

## GET A PART EXCHANGE VALUATION

Enter your details below to get a valuation for your vehicle in part exchange.

## LEARN MORE ABOUT THIS VEHICLE

FREEPHONE:

(+44) 0330 1747 555

brotherwood.com



Contact us today for more information and to discuss your needs. Our dedicated, experienced team are on hand to answer all of your questions about Wheelchair Accessible Vehicles.

★★★★★ 315 reviews