



MERCEDES-BENZ V300 D EXCLUSIVE LONG KLASTAR [PO4285]

Price on Application



DETAILS

Headroom (in.)	57.1
Transmission	<u>Automatic</u>
Mileage	7399
Seats	4 + Wheelchair
Model Year	2025
Colour	Obsidian Black
Engine	2.0D
Fuel Type	Diesel
Brand	Mercedes-Benz

The Brotherwood® 'Klaster' Mercedes-Benz V-Class Exclusive WAV conversion is designed and engineered to be the ultimate Wheelchair Accessible Vehicle. This 2025 model is highly equipped with Parking Package, 360 Degree Camera, Driving Assistance, Memory Pack, Burmester Surround Sound, Fridge Pack and Smartphone Integration.

SKU: PO4285

We're proud
to be an
RAC Approved
Dealer

RAC
Approved Dealer



VEHICLE DESCRIPTION

ULEZ Compliant

Vehicle Summary

Dimensions

	mm	inches
Entry Height	1425-1474	56.1-58.1
Max. Internal Headroom	1450	57.1
Max. Wheelchair Armrest Width	660	26
Lowered Floor Length	2200	86.5
Lowered Floor Width	711	28
Access Ramp Length	1070	42
Access Ramp Width	740	29
Access Ramp Angle		12.7°
Exterior Length	5140	202.4
Exterior Width with mirrors	2249	88.5
Exterior Height	1880	74

Internal headroom may vary depending on wheelchair travelling position.

Every effort has been made to ensure the accuracy of all information. However, errors may occur. Please check with the sales team about items which may affect your decision to purchase. The majority of our customers are exempt VAT, hence all prices shown exclude VAT.



Conversion Features

- Powered Access Ramp
- High-Traction Ramp Surface
- 2 x Second-Row Seat
- LED Lighting Kit in Rear
- Sculpted Headlining for +2" Entry Height
- Mercedes-Benz Illuminated Stainless Steel Rubbing Strips
- Four Point Webbing Restraint
- Quattro Rear Restraints
- Lowered, Level Floor
- Hook-i PAWRS (combined winch & front restraints)

Specification Includes

- Parking Package with 360 Degree Camera
- Attention Assist
- Cruise Control
- Smartphone Integration Package
- Driver Assistance Pack
- Memory Pack
- Burmester Surround Sound System
- Refrigerator Pack



- 19" AMG alloy wheels
- Automatic powered tailgate with separate opening rear window
- LED Intelligent Light System
- Heated windscreen wash
- Privacy glass
- Ambient lighting - 64 colours
- Speakers with Fontbass
- Electric sliding doors - both sides
- Heated, lumbar support comfort front seats, including armrests
- Luxury centre console with retractable cover, two cup holders, 12 V socket and multimedia box
- Silk Beige Nappa leather upholstery
- MBUX Infotainment system with voice control for navigation, phone calls, music and climate adjustment

GET A PART EXCHANGE VALUATION

Enter your details below to get a valuation for your vehicle in part exchange.

ASK US A QUESTION

FREEPHONE:

[\(+44\) 0330 1747 555](tel:+441747555555)

Contact us today for more information and to discuss your needs. Our dedicated, experienced team are on hand to answer all of your questions about Wheelchair Accessible Vehicles.